

## **Product Information**

## Protein Endotoxin Removal Kit

## CAT#: EKIT-040

Protein Endotoxin Removal Kit, also known as protein sample endotoxin removal kit or Protein Endotoxin Removal Kit, is a kit that can efficiently remove more than 99% of the endotoxin in protein solution samples within 15 minutes. Protein samples with an initial endotoxin concentration of 100EU/ml can be treated with this kit with endotoxin levels below 0.1EU/ml to a level considered free of endotoxin contamination. Therefore, this kit is widely used in endotoxin treatment of samples involving cell or animal experiments.

## **Product Information**

similar	EKIT
Alias	EKIT
productkeywords	Protein Endotoxin Removal Kit; EKIT-040
Applications	This kit is optimized for the removal of endotoxin from protein samples and can also be used for the removal of endotoxin from nucleic acid or small molecule solution samples, but the removal of endotoxin from nucleic acid samples is slightly less effective than that of protein samples.
Storage	Store at -20°C for two years. Store at 4°C for 3 months.