

Product Information

Endotoxin Indicator

CAT#: E027

The endotoxin indicator is used for GMP verification of dry heat sterilization equipment and endotoxin interference test of medical equipment in pharmaceutical production enterprises, confirming the heat distribution at different positions in tunnel oven and dry heat oven chamber, confirming whether the original program of the scheduled sterilization and heat removal can meet the pre-design requirements.

CAT#	PRODUCT	Package	Storage
E027	Endotoxin Indicator	2 mg×20	2-8 °C

Product Information

Forms	Lyophilized powder
Storage	2-8 °C
Principle	Gel Clot LAL
Origin	E.coli O111:B4
Citations	Use endotoxin free water to configure